

Bring Haemostasis to light



## Breakthrough innovation in Thrombin Generation (TG) ST Genesla: A complete solution

100% standardised

100% automated

- **Specific features:** 6 TG parameters + ETP inhibition
  - Both in absolute & normalised units
  - · All parameters calculated automatically.
- Assayed Reference Plasma for result normalisation.
- 3 QC levels to cover the entire working range.
- **Trigger reagents and QC combined** for improvement of standardisation.
- Precise temperature control at 37°C.

- Routine features: Calibration & QC management (incl Westgard), complete traceability and STAT samples.
- User friendly interface with secured access.
- Positive identification & continuous loading of reagents, samples & disposables.
- Unitary cuvettes preloaded on trays.

Walk-away benchtop analyser

100% innovative

- New patented calibration method: only once daily required, insensitive to anticoagulant drugs.
- Optimised reagents for hypo-, hyper-coagulable & samples containing anticoagulant.
- Protein C pathway function assessed by Thrombomodulin reagent.
- Ready to use fluorescent substrate.



Diagnostica Stago S.A.S. RCS Nanterre B305 151 409 3, allée Thérésa 92600 Asnières sur Seine France Ph.: +33 (0)1 46 88 20 20 Fax: +33 (0)1 47 91 08 91

Fax: +33 (0)1 47 91 08 91 webmaster@stago.com